

Select™
GC Columns

Select Al₂O₃ MAPD : a new Inert Al₂O₃
column for the trace analysis of polar
reactive hydrocarbons.

Jaap de Zeeuw and Coen Duvekot
Varian, Inc.
Middelburg, The Netherlands



Varian

PO Box 8033, 4330 EA Middelburg
The Netherlands

Phone : 31-118-671000

Fax : 31-118-671325

Email : North and South America

Europe, Asia, Africa & Australia

: helpdesk.us@varianinc.com

: helpdesk.eu@varianinc.com



P-154

Select Al₂O₃ MAPD : a new Al₂O₃ column for the trace analysis of polar reactive hydrocarbons.

Jaap de Zeeuw and Coen Duvekot
Varian, Inc. Middelburg, The Netherlands

Introduction

In petrochemical arena analysis of volatile hydrocarbons and impurities has become a routine application. In the C₁-C₅ range many unsaturated hydrocarbons are present (alkenes, alkynes) and also these have to be separated and quantified. As hydrocarbons are generally believed to be inert they will elute from a wide variety of stationary phases. Liquid type phases like 100% polydimethyl siloxanes will separate the saturated hydrocarbons, but selectivity is not present for separating the unsaturated hydrocarbons. Highly selective materials have been developed for the separation of unsaturated hydrocarbons and especially the adsorption chromatography has proven to be very effective. Materials like alumina [1-3], silica[4] and carbon[5] have unique selectivities and are widely used. Also for the separation and quantification of polar impurities in hydrocarbons streams adsorbents like PoraBOND, PoraPLOT and Lowox have found wide application. Most widely used is the aluminium oxide. Aluminium oxides separate all C₁-C₅ hydrocarbons but has to be deactivated. The current alumina column show a lack of inertness towards reactive hydrocarbons *cq*, propadiene and methyl acetylene. This translates in a reduced response and a non reproducible response in time.

A new Al₂O₃ adsorbent is introduced called Select Al₂O₃ MAPD. The new material is intensively deactivated which provides highest response for reactive hydrocarbons and also generates a stable response in time.

Aluminium oxide

Alumina adsorbents in capillary columns have been introduced in 1963 and were commercialized in fused silica capillary columns in 1981. The alumina adsorbent has a very high activity and will retain components as light as ethane. To make the highly active alumina work as stationary phase in gas chromatography, it has to be deactivated. Deactivation can be done in many ways; however, the most practical and reproducible way is the deactivation with inorganic salts. Very popular is the KCl deactivated Alumina which provides a general non-polar alumina surface. Such a column will elute acetylene before the butane as seen in Figure 1. The alumina surface can be made more polar by deactivation with a salt

with a higher charge, for instance sodium sulfate. The resulting alumina layer will elute the acetylene after the butane peaks indicating the higher polarity., see fig.2. Also the separation of trans-2-butene and 1-butene is much improved which is especially useful for impurity analysis in butene streams.

The selectivity of alumina for hydrocarbons is very high. All C₁-C₄ hydrocarbons can be separated baseline. The resolution between the different hydrocarbons is sufficient to measure many traces of C₁-C₄ hydrocarbons in a matrix of any C₁-C₄ hydrocarbon. For this reason alumina is one of the most widely used columns in petrochemistry for analyzing hydrocarbon impurities in various light hydrocarbon matrices.

Alumina columns have been developed in all dimensions, from 0.25 mm ID up to 0.53 and also in UltiMetal tubing.




Figure 1: C₁ - C₆ hydrocarbons on Al₂O₃, KCl deactivated adsorbent PLOT
 Column: 50 m x 0.32 mm fused silica, Al₂O₃/KCl df=5 um; Oven: 70 °C(4 min)-->180 °C, 10°C/min;
 Carrier gas: Helium; Injection:Split; Detection: FID; 1:methane; 2:ethane; 3:ethylene; 4:propane;
 5:cyclopropane; 6:propylene; 7:acetylene; 8:isobutane; 9:propadiene; 10:butane; 11:trans-2-butene; 12:1-butene; 13:isobutene; 14:cis-2-butene; 15:isopentane; 16:methylacetylene; 17:pentane; 18: 1,3-butadiene; 19:ethylacetylene; 20:hexane




Figure 2: C₁ - C₄ hydrocarbons on Al₂O₃, Na₂SO₄ deactivated adsorbent PLOT
 Column: 50 m x 0.32 mm fused silica, Al₂O₃/KCl df=5 um; Oven: 120 °C; Carrier gas: Helium;
 Injection:Split; Detection: FID; 1. Methane; 2. Ethane; 3.Ethene; 4.Propane; 5.Propene;
 6.Butane; 7.Propadiene; 8.Acetylene 9. tr-2-Butene;; 10. Butene-1; 11. iso-Butene; 12. cis-2-Butene; 13. 1,3-Butadiene; 14.Methylacetylene;

Problems with alumina

Although alumina has unique separation characteristics, it also has limitations. The activity of the adsorbent will adsorb any moisture, carbon dioxide or other polar impurities in the sample. If moisture is brought on alumina, the retention

times for hydrocarbons start to decrease. The reduction of retention time is due to water deactivating the sorbent involved. Water can simply be removed by heating up the column to 200 °C for 10-15 minutes. Water will elute and the column is then regenerated. If isothermal set-up is required, one could use a polar pre-column to retain water. A polyethylene glycol coated column with a 1.2 micron film works very well as the C1-C6 hydrocarbons will elute from this column before water elutes. Compounds with some degree of polarity (like a functional group) do not elute or are strongly retained (alcohols, esters ethers, bases, acids). Alumina has been used for analysis of CFC (Chloro-Fluoro-Carbons) which only worked for the fully substituted molecules. When a hydrogen

was present the molecule was adsorbed.




Fig. 3

Practically the application of Alumina is restricted to hydrocarbons.

Also on hydrocarbon analysis several unwanted interactions can occur. Figure 3 shows a reactivity issue when analyzing propadiene and methyl acetylene.




Fig. 4

Figure 4 shows the behavior of an aluminium oxide column when used in a process environment running then same analysis over a period of 6 months. The response of propadiene relative to propane is plotted. During this period the response decreased with near 30%.




Fig. 5

The absolute response to the reactive hydrocarbons varied strong from column to column and was much lower then theoretically would be expected. The

typical response for all hydrocarbons is shown in figure 5: Propadiene (PD) response is approx. 35% and methyl acetylene (MA) is approx. 70%. Table 1 show the response factor for MA and PD relative to propane of 10 commercially available Al₂O₃ columns.

Table 1 : Response of MA and PD relative to propane

	Response MA (%)	Response PD (%)
Column 1 supplier 1	59	24
Column 2 supplier 1	63	24
Column 3 supplier 1	59	18
Column 1 supplier 2	73	66
Column 2 supplier 2	73	57
Column 1 supplier 3	64	32
Column 2 supplier 3	56	27
Column 3 supplier 3	63	24
Column 4 supplier 3	59	25
Select Al ₂ O ₃ MAPD	93	84

The reactivity towards these hydrocarbons is related to their elution temperature and retention time. This is not a surprise, as longer elution times gives the component more chance to interact and also higher elution temperatures will not contribute to a better response.

Figure 6 shows the response of different hydrocarbons run on a polar type Al₂O₃ capillary but the analysis was done using different column flow rates. With increased flow rate the response also increases significantly for propadiene and methylacetylene.




Fig. 6

A new Al₂O₃ column: Select Al₂O₃ MAPD

As it is not possible to operate columns at such high flow rates there is a clear room for improvement on the Al₂O₃ and challenge was to improve the response and stability for the reactive hydrocarbons methyl acetylene and propadiene.

The new alumina was given a special treatment and was deactivated much more intense than usually done with Al₂O₃/KCl and the Al₂O₃/Na₂SO₄ adsorbents.

As a result the retention and reactivity was decreased and a highly stable column was obtained especially for methyl acetylene(MA) and propadiene(PD). The selectivity of the Select Al₂O₃ MAPD is close to the Na₂SO₄ selectivity

although it is not exactly

identical and small variations can occur. The response for propadiene is more than doubled, see figure 7.

Typical characteristics of this new column are:




Fig. 7

Up to a factor 2 higher response for acetylene, methylacetylene and propadiene;

This is especially important when running traces. In many industrial hydrocarbon streams there are tight specifications on these reactive




Fig. 8

hydrocarbons. Specification which are in the order of 1 – 3 ppm in a ethylene, propylene, butylene and butadiene product. Figure 8 shows an analysis of a 1,3 Butadiene stream. Note the high responses for the reactive hydrocarbons.

Improve detection limits up to a factor 2-3

As the signal produced for a propadiene peak is 2-3 times higher, so is the signal – noise ratio. As a result detection limits can be lowered with a similar factor. This makes the method more accurate and will also provides more accurate measurements.

No change of response in time

The Select Al₂O₃ MAPD was tested in a on-line analyzer measuring continuously light hydrocarbons including MA and PD over a period of 18 weeks, see fig.9 No change in response factor was observed during this period which was a welcome improvement for this hydrocarbon method.




Fig. 9

Shorter run times

The run times for most hydrocarbon applications are shorter due to the lower retention of the column. Although the amount of Aluminium oxide in the column is the same, only the least active sites will participate in the retention process(the sites with highest activity are all deactivated). As a result the peak shape improves and peaks elute faster. The loadability was not found to be significantly lower.

Impact of water

Despite of the better inertness the Select Al₂O₃ MAPD remains sensitive to water in the sample or the carrier gas. Water will be trapped on the alumina and will make components elute faster. Regeneration(water bake-out) is done by baking 5-10 minutes at 200 C. This can be done periodically or after each analysis. Regeneration will be approximate 2 times faster with the Select Al₂O₃ MAPD capillary. If an isothermal setup

is preferred one can take a CP-Wax thick film pre-column, let the c1-C6 hydrocarbons go through and retain the water on the CP-wax column. After hydrocarbons are lead into the alumina column, a valve switch sent the later eluting water peak to vent or to the detector.

Application

This new alumina column can be used for all applications where C1-C10 hydrocarbons are analyzed and where there is an interest in quantitative analysis of reactive hydrocarbons like propadiene and methyl acetylene. Also compounds like acetylene, ethyl - and vinylacetylenes have least chance for reactions using this inert alumina.

We have seen good performance for % and impurity analysis see figures 10-14.




Fig.10: Impurities in propylene; 1 Methane; 2 Ethane; 3 Ethylene; 4 Propane; 5 Propylene; 6 Cyclopropane; 7 Butane 8 Propadiene; 9 Acetylene; 10 Trans-2-Butylene; 11 Butylene; 12 Iso-Butylene; 13 Cis-2-Butylene; 14 1,3-Butadiene; 15 Methylacetylene;




Fig 11 Elution of acetylenes

Column : 50 m x 0.53 mm Select Al₂O₃ MAPD; Oven: 40°C, 5min → 160°C, 10°C/min → 200°C, 20°C/min; Carrier gas : 4 psig, 4 min → 11 psig, 0.5 psig/min

1 Ethane; 2 Propane; 3 Trans-2-butene; 4 Iso-pentane; 5 Methyl acetylene; 6 Vinyl acetylene 7 Ethyl acetylene; Ref: Jim Luong, Dow Chemical Canada




Fig.12
 Column : 50 m x 0.53 mm Select Al₂O₃ MAPD; Oven: 40°C, 5min → 160°C, 10°C/min → 200°C, 20°C/min; Carrier gas : 4 psig, 4 min → 11 psig, 0.5 psig/min
 1 Methane; 2 Ethylene; 3 Iso-butane; 4 n-butane; 5 Iso-pentane; 6 Cis-2-butene; 7 n-pentane
 ref: Jim Luong, Dow Chemical Canada




Fig.13
 Column : 50 m x 0.53 mm Select Al₂O₃ MAPD; Oven: 40°C, 5min → 160°C, 10°C/min → 200°C, 20°C/min; Carrier gas : 4 psig, 4 min → 11 psig, 0.5 psig/min
 1 Propylene; 2 Propadiene; 3 Trans-2-butene; 4 Iso-butene; 5 Cis-2-butene; 6 n-pentane; 7 1,2 butadiene; 8 1,3-butadiene
 ref: Jim Luong, Dow Chemical Canada




Fig.14

Column : 50 m x 0.53 mm Select Al₂O₃ MAPD; Oven: 40°C, 5min → 160°C, 10°C/min → 200°C, 20°C/min; Carrier gas : 4 psig, 4 min → 11 psig, 0.5 psig/min

1 Methane; 2 Ethane; 3 Ethylene; 4 N-butane; 5 Propadiene; 6 1-butene; 7 Iso-butene
8 1,2-butadiene; 9 1,3-butadiene; 10 Ethyl acetylene

ref: Jim Luong, Dow Chemical Canada

Conclusion

A new Al₂O₃ PLOT columns is introduced with a better deactivation. The Varian Select Al₂O₃ MAPD shows improved inertness towards reactive hydrocarbons, especially propadiene and methylacetylene when analyzed down to ppm levels. Response for these components is much higher and also the response factor in time remains constant making the column very reliable in a analyzer environment.. The Select Al₂O₃ MAPD separates light hydrocarbons in the C1-C10 range can be used for similar applications as the existing Al₂O₃ columns. Nett benefit will be better detectivity and a consistent data.

References

- 1 .J.Kirkland, Anal. Chem., 35(1963)1295-1297;
- 2 W.Schneider, J.C.Frohe, H.Brudereck, J.Chromatogr., 155(1978)311-327;
- 3 R.de Nijs and J.de Zeeuw, J.Chromatogr., 279(1983)41-48;
- 4 M..Mohnke, D.Estel, J.de Zeeuw and C.Duvekot, American Laboratory, may 1997, p.22-24;
- 5 J. de Zeeuw,Encyclopedia of Separation Science,Gas Solid Chromatography, 1992, p.481-488;

Acknowledgements

The authors want to acknowledge the support from Jim Luong and Ronda Gras, Dow Chemical Canada, for generating test data and stimulating discussions.

Ordering Information

Capillary columns

Varian Select Al₂O₃ MAPD 50 m x 0.32 mm	CP7431
Varian Select Al₂O₃ MAPD 50 m x 0.53 mm	CP7432

Fused silica coupling

Universal CP-QuickSeal	(10pc)	CP4787
Universal CP-QuickSeal Splitter	(5pc)	CP4798

Reamers

For graphite and vespel ferrules	CP8427
----------------------------------	--------

Ferrules

Graphite ferrules 1/16" for 0.53 mm	(10pc)	CR211108
Graphite ferrules 1/16" for 0.53 mm in HP-series	(10pc)	CR211168
Graphite ferrules 1/16" for 0.32 mm	(10pc)	CR211105
Graphite ferrules 1/16" for 0.32 mm in HP-series	(10pc)	CR211165
Graphite Vespel ferrules 1/16" for 0.32 mm	(10pc)	CR213105
Graphite Vespel ferrules 1/16" for 0.32 mm in HP series	(10pc)	CR213165
Graphite Vespel ferrules 1/16" for 0.53 mm	(10pc)	CR213108
Graphite Vespel ferrules 1/16" for 0.53 mm in HP series	(10pc)	CR213168

Gas Leak detector

220 V Eur	CP80610
240 V UK	CP83610
110 V USA	CP87610

Filters

Gas-Clean Oxygen filter	CP17970
Gas-Clean Moisture filter	CP17971
Gas-Clean Charcoal filter	CP17972
GC-MS Filter	CP17973
Connecting unit for 1 filter 1/4" tubing	CP7980
Connecting unit for 1 filter 1/8" tubing	CP7988
Connecting unit for 2 filter 1/4" tubing	CP738406
Connecting unit for 2 filter 1/8" tubing	CP738407
Connecting unit for 4 filter 1/4" tubing	CP7989
Connecting unit for 4 filter 1/8" tubing	CP736520
Flush head for Connection unit	CP7987

Column Cutting

Tungsten Carbide Glass Knife	CP22674
Carborundum Glass file	CP22672

Varian, Inc. is a high-tech organisation with head quarter in Palo Alto, California, USA. It is a major supplier of instruments and products for GC, GC/MS, HPLC, SPE, AA, ICP/MS, UV-VIS, NMR and vacuum products.